

Preanalytical Considerations about Circulating Tumor DNA

Specimen Improvements, Standards & Technologies



LIQUID BIOPSY SUMMIT

20 – 21 February 2020 Lisbon, Portugal

Dr. Daniel Grölz

- Need for pre-analytical workflow standardization to reduce diagnostic errors
- Pre-analytical factors that influence the outcome of ccfDNA analysis
- International initiatives and requirements to standardize pre-analytical workflows
- Liquid biopsy preservation and workflow solutions

- Need for pre-analytical workflow standardization to reduce diagnostic errors
- Pre-analytical factors that influence the outcome of ccfDNA analysis
- International initiatives and requirements to standardize pre-analytical workflows
- Liquid biopsy preservation and workflow solutions



- ❖ 150,000 papers documenting thousands of claimed biomarkers, but fewer than 100 have been validated for routine clinical practice

Bring on the biomarkers, George Poste, Nature 2011

- ❖ Diagnostic errors cause about 10% of all patient deaths and about 17% of adverse events

Institute of Medicine (IOM) Report Sept. 2015

- ❖ The pre-analytical phase accounts for 46% to 68% of errors observed during the total testing process

Medical Laboratory Observer, May 2014

- ❖ Unnecessary expenditure caused by pre-analytical errors in a typical U.S. hospital (~ 650 beds) of ~ \$1.2 million per year

Green SF. Clin Biochem. 2013

- ❖ Irreproducible preclinical research exceeds 50%, US \$28B / year spent on preclinical research that is not reproducible - in the US

Freedman LP, Cockburn IM, Simcoe TS (2015) PLoS Biol 13(6): e1002165.doi:10.1371/journal.pbio.1002165

An Analytical Test Result is the Result of an Entire Workflow




Specifying, developing and verifying preanalytical workflows has to be part of the analytical test development



European Conference. Standards: Your Innovation Bridge. Brussels (2014). SPIDIA Booth.


Provided by courtesy of SPIDIA

Transcription induction and degradation in human EDTA Blood stored at Room Temperature



IL-1b mRNA

Guenther K. et al. AMP Poster (2005)



c-fos mRNA

Guenther K. et al. CLI 5, 26-28 (2008)

Provided by courtesy of D. Groelz, PreAnalytiX

- Need for pre-analytical workflow standardization to reduce diagnostic errors
- **Pre-analytical factors that influence the outcome of ccfDNA analysis**
- International initiatives and requirements to standardize pre-analytical workflows
- Liquid biopsy preservation and workflow solutions


Pre-analytical Factors that Influence the Outcome of ccfDNA Analysis

Preanalytical step affecting ccfDNA	Challenge	Recommendation
Blood collection tube	Release of genomic DNA from leukocytes; PCR inhibition	Use of a dedicated ccfDNA stabilization tube Use of EDTA tube, in case no dedicated ccfDNA stabilization tube is available
Time between blood collection and plasma processing	Release of genomic DNA from leukocytes	Must be validated in combination with downstream application In case of unstabilized EDTA blood within 2 to 6 hours
Plasma or serum	Release of genomic DNA from leukocytes	Use of plasma
Plasma processing protocol	Incomplete separation of cellular fraction Mechanical lysis of blood cells	For EDTA blood use double centrifugation protocol with low and high speed centrifugation Follow recommendations of ccfDNA tube manufacturer
Plasma storage	Reduction of yield Increased fragmentation	Do not store plasma at 4°C for longer than 24 hours Freeze at -20°C or -80°C Avoid repeated freezing/ thawing cycles
DNA purification method	Suboptimal compatibility with blood collection tube Low yield	Validation of DNA purification method Follow recommendations of ccfDNA tube manufacturer
DNA quantification	Over-quantification due to absorption of impurities in spectrophotometry	Validation of QC-methods; use of qPCR based methods
DNA storage	Reduction of yield Increased fragmentation	Store ccfDNA at -20°C or below Avoid repeated freezing and thawing

Grölz et al. *Curr Pathobiol Rep.* 2018;6(4):275-286

Preanalytical step affecting ccfDNA	Challenge	Recommendation
Blood collection tube	Release of genomic DNA from leukocytes; PCR inhibition	Use of a dedicated ccfDNA stabilization tube Use of EDTA tube, in case no dedicated ccfDNA stabilization tube is available
Time between blood collection and plasma processing	Release of genomic DNA from leukocytes	Must be validated in combination with downstream application In case of unstabilized EDTA blood within 2 to 6 hours
Plasma or serum	Release of genomic DNA from leukocytes	Use of plasma
Plasma	<p>„the entire workflow, including specimen/sample storage and transport conditions, and its impact on the stability of biomolecules intended to be examined shall be verified and validated“</p> <p>ISO 20186-3:2019 Molecular in vitro diagnostic examinations — Specifications for pre- examination processes for venous whole blood — Part 3: Isolated circulating cell free DNA from plasma.</p>	
Plasma		
DNA		
	Low yield	Follow recommendations of ccfDNA tube manufacturer
DNA quantification	Over-quantification due to absorption of impurities in spectrophotometry	Validation of QC-methods; use of qPCR based methods
DNA storage	Reduction of yield Increased fragmentation	Store ccfDNA at -20°C or below Avoid repeated freezing and thawing

Apoptosis of white blood cells leads to release of high molecular weight DNA




ccfDNA extracted from plasma collected using EDTA tubes from a subject directly after blood draw (t0) and after seven days at room temperature (T7d, 25°C)

One-microliter eluates analyzed using Agilent High Sensitivity DNA Kit

Provided by courtesy of D. Groelz, PreAnalytiX

Apoptosis of white blood cells leads to release of high molecular weight DNA




ccfDNA extracted from plasma collected using EDTA tubes from a subject directly after blood draw (t0) and after seven days at room temperature (T7d, 25°C)

One-microliter eluates analyzed using Agilent High Sensitivity DNA Kit

Provided by courtesy of D. Groelz, PreAnalytiX

Apoptosis of white blood cells leads to increased yield and dilution of the target ccfDNA



X whole blood storage duration at room temperature before plasma generation (t1,3,6: 1, 3, and 6 days)

Y ratio of 18S rDNA copy numbers determined in plasma after indicated blood storage durations versus immediately after blood collection (t0)


ISO 20186-3:2019

Molecular in vitro diagnostic examinations — Specifications for pre-examination processes for venous whole blood
— Part 3: Isolated circulating cell free DNA from plasma. Annex A.

Follow manufacturer's centrifugation protocols to avoid gDNA carryover or reduced volume of separated plasma


Manufacturer recommendation for PAXgene Blood ccfDNA Tube

- 1st Spin 15 min at 1600-3000 x g
- 2nd Spin 10 min at 1600-3000 x g
- Medium brake




Best outcome: Maximum recovery of clear plasma with compact level interface


Reduced
centrifugation time:
5 min 1600 x g



Maximum
brake:




Provided by courtesy of D. Groelz, PreAnalytiX



Blood from 8 donors collected in EDTA and PAXgene Blood ccfDNA Tubes, manual (QIAamp cNA kit) and automated ccfDNA (QIAasymphony) isolation;
DNA quantification with Cell-free DNA ScreenTape Assay (4200 TapeStation System), High Sensitivity DNA Analysis Assay (Bioanalyzer 2100), Qubit® 1x (premixed) and 10x (dilution needed) dsDNA HS Assay (Qubit 2.0), Investigator Quantiplex® Pro Assay (QIAGEN Rotor-Gene® Q), 18S rDNA 66 bp (inhouse assay on QX200 Droplet Digital™ PCR).

Provided by courtesy of
D. Groelz, PreAnalytiX

Post Blood Collection ccfDNA Profile Changes - Impact on EGFR Test



X venous whole blood storage duration (in days) before plasma preparation

Y $\Delta CT = CT (\text{mutant}) - CT (\text{wildtype control})$

■ EDTA Blood

■ Stabilized Blood

.... Threshold (given by the examination provider)

The average of 7 donors is shown

ISO 20186-3:2019

Molecular in vitro diagnostic examinations — Specifications for pre-examination processes for venous whole blood

— Part 3: Isolated circulating cell free DNA from plasma. Annex A.

Blood collected in EDTA and PAXgene Blood ccfDNA tubes

Spiked with restriction enzyme treated EGFR DNA with mutation T790M, equivalent to 200 copies

ccfDNA tested with the commercially available EGFR Plasma PCR Kit (RUO)

- Need for pre-analytical workflow standardization to reduce diagnostic errors
- Pre-analytical factors that influence the outcome of ccfDNA analysis
- International initiatives and requirements to standardize pre-analytical workflows
- Liquid biopsy preservation and workflow solutions

SPIDIA – FP7 (2008 – 2013)



- ⇒ 16 Partners
- New technologies for sample collection, stabilization, processing, transport, storage (Blood, Tissues)
- 9 EU CEN Standards

SPIDIA4P – H2020 (2017 – 2020)



- ⇒ 19 Partners
- ⇒ 14 associated consortia & stakeholder organizations
- 13 additional new CEN & ISO Standards
- EQAs
- European implementation

www.spidia.eu ⇒ **New Website. Subscribe the Newsletter!**

The SPIDIA project has received funding under the Seventh Research Framework Program of the European Union, FP7-HEALTH-2007-1.2.5, under grant agreement no. 222916. The SPIDIA4P project receives funding from the European Union's Horizon 2020 research and innovation programme under grant agreement no. 733112.



**Vienna
Agreement
1991**



- 2019: 8 ISO/International Standards
- 2014: 8 new projects for ISO Standards approved in ISO/TC 212 „Clinical laboratory testing and in vitro diagnostic test systems”



- 2015: 9 CEN Technical Specifications published
- 2013: 9 new projects approved in CEN/TC 140 „In vitro diagnostic medical devices“
- 2010: Start of standardization work



*European Conference. Standards:
Your Innovation Bridge. Brussels
(2014). SPIDIA Booth.*





1. Problem - Errors in Diagnostics





■ Preanalytics ■ Analytics ■ Postanalytics

2. Technical Solutions



3. Ring-Trials – Blood RNA (l.) and DNA (r.)






Provided by courtesy of SPIDIA

❖ CEN

- ❖ Recognized by the EU and the European Free Trade Association (EFTA) as being responsible for developing standards at European level
- ❖ Development of a European Standard (EN) or International Standard (ISO) is governed by the principles of consensus, openness, transparency, national commitment and technical coherence
- ❖ One European Standard replaces 34 national standards


❖ CEN/TC 140 (Committee for in vitro diagnostic medical devices)

- ❖ 34 EU countries National Standards Bodies (NSB)
- ❖ Stakeholders in liaison & cooperations
 - ❖ **European Commission** (EC), **ESP** (European Society of Pathology), **EFLM** (European Federation of Laboratory Medicine), **IFCC** (Int. Federation of Clinical Chemistry and Laboratory Medicine), **JISC** (Japanese Industrial Standards Committee), **MedTech** (Alliance of European medical technology industry associations, founded by **EDMA**), **EPBS** (European Association for Professions in Biomedical Science), **BBMRI-ERIC** (Biobanking and BioMolecular resources Research Infrastructure - European Research Infrastructure Consortium), in progress, **ISO/TC 212** (Clinical laboratory testing and in vitro diagnostic test systems), **ISO/TC 276** Biotechnology



● Molecular in-vitro diagnostic examinations - Specifications for pre-examination processes for

- Blood — Cellular RNA
- Blood — Genomic DNA
- Blood — Circulating cell free DNA
- FFPE tissue — DNA
- FFPE tissue — RNA
- FFPE tissue — Proteins
- Frozen tissue — RNA
- Frozen tissue — Proteins
- Metabolomics in urine, serum and plasma



Provided by courtesy of SPIDIA

Molecular in-vitro diagnostic examinations - Specifications for pre-examination processes for

- Blood — Cellular RNA, gDNA, ccfDNA, ccfRNA
- Blood – *Exosomes, ccfRNA*
- Circulating Tumor Cells – **DNA, RNA, staining**
- Tissue (FFPE) — DNA, RNA, Proteins
- Tissue (Frozen) – RNA, Proteins, DNA
- Tissue (FFPE) - *staining*
- Fine Needle Aspirates – *DNA, RNA, Proteins*
- Saliva – **DNA**
- Urine & Body Fluids – *cfDNA*
- Microbiome – *Stool, Saliva etc.*

published CEN
published ISO
in development


<https://www.spidia.eu/>




- New European In Vitro Diagnostic Regulation in force since May 2017
- Pre-analytical workflow parameters in several sections
 - 6. PRODUCT VERIFICATION AND VALIDATION (Annex II)
 - 6.1. Information on analytical performance of the device
 - 6.1.1. Specimen type

This Section shall describe the different specimen types that can be analysed, including their stability such as storage, where applicable specimen transport conditions and, with a view to time-critical analysis methods, information on the timeframe between taking the specimen and its analysis and storage conditions such as duration, temperature limits and freeze/thaw cycles

PreAnalytiX a QIAGEN/BD Company and SPIDIA4P partner



**Product
Development
Process**



Certification according
to ISO 13485

Company Quality Manual:
Process Landscape

Global Process
SOPs incl. legal
requirements

Technical SOPs for pre-analytical workflows
based on ISO & CEN standards

Provided by courtesy of D. Groelz, PreAnalytiX

Contents	Page
Foreword	iv
Introduction	v
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 General consideration	5
5 Outside the laboratory	5
5.1 Specimen collection	5
5.1.1 Information about the specimen donor/patient	5
5.1.2 Selection of the venous whole blood collection tube by the laboratory	6
5.1.3 Venous whole blood collection from the donor/patient and stabilization procedures	6
5.1.4 Information about the specimen and storage requirements at the blood collection facility	7
5.2 Transport requirements	7
6 Inside the laboratory	8
6.1 Specimen reception	8
6.2 Storage requirements for blood specimens	8
6.3 Plasma preparation	9
6.4 Storage requirements for plasma samples	9
6.5 Isolation of the ccfDNA	10
6.5.1 General	10
6.5.2 Using blood collection tubes with stabilizers	10
6.5.3 Using blood collection tubes without stabilizers	11
6.6 Quantity and quality assessment of isolated ccfDNA	11
6.7 Storage of isolated ccfDNA	11
6.7.1 General	11
6.7.2 ccfDNA isolated with commercially available kits	12
6.7.3 ccfDNA isolated with the laboratory's own protocols	12
Annex A (informative) Impact of pre-examination process steps on circulating cell free DNA profiles in venous whole blood plasma	13
Bibliography	16

ISO 20186-3:2019 - Molecular in vitro diagnostic examinations — Specifications for pre-examination processes for venous whole blood — Part 3: Isolated circulating cell free DNA from plasma


- Need for pre-analytical workflow standardization to reduce diagnostic errors
- Pre-analytical factors that influence the outcome of ccfDNA analysis
- International initiatives and requirements to standardize pre-analytical workflows
- Liquid biopsy preservation and workflow solutions

Pre-analytical Steps: Part of a Whole Diagnostic Test Workflow



Provided by courtesy of D. Groelz, PreAnalytiX

Multimodal Analysis: Extraction and Analysis of CTC RNA, ccfDNA and gDNA from a Single Blood Sample



*Babayan et al. [abstract].
Proceedings: AACR Advances in
Liquid Biopsies 2020;
January 13-16, 2020; Miami, FL.*

* For Research Use Only. Not for use in diagnostic procedures.

- Unless otherwise stated, all products mentioned in this presentation are for Research Use Only and not for use in diagnostic procedures.
- **RUO** products shown here are for Research Use Only. Not for use in diagnostic procedures.
- **MBA** products shown here are for Molecular Biology Applications Only. Not for the diagnosis, prevention or treatment of a disease.
- For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN or PreAnalytiX kit handbook or user manual. QIAGEN or PreAnalytiX kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.
- A CE-marked PAXgene Blood ccfDNA Tube for in vitro diagnostics is available in certain European countries. The performance characteristics of the RUO and the IVD tubes are equivalent.